

# homato

## HOME MAGNETOTHERAPY SYSTEM



MAINS

0 20 40 60 80 100

B

0.1  
KL 2.5



START  
STOP

1

Norra Medicin

# [H]O[M]A[G]N[E]T[O]THERAPY SYSTEM

Homato is the acronym of **H**ome **M**agnet**o**therapy System, an analogue model in the Norra Medicin family of products. It is designed to be used by non medical-doctor practitioners. This is a high-intensity magnetic induction field system, an effective therapy for the treatment of musculoskeletal pain and inflammation.

Homato has a modular but non-algorithmic program, running for 1-30 minutes, depending.

Homato has an adjustable intensity, with levels from 2 to 100.


At the lowest intensity, the pulse rate is close to 2 pulses per second and at the highest intensity level the rate can go up to 300 per second.

## Parameters


Oscillation frequency	2 to 50 Hz (S, 50kHz)
Oscillation	> 10
Magnetic induction	1000 mT
Performance (dB max)	100.500 T/s
Program Duration	10-20 minutes
Input Power	230 V AC
Unit Weight	1+5 kg
Unit Dimensions	38cm x 32cm x 10.5cm
Applicator diameter	
Warranty, control unit	4 years
Warranty, applicator	5 yearss
Return Period	30 days
Restocking Fee	20% of purchase price

## Benefits


- Easy to use, high comfort product.
- Five modulation options for different conditions.
- No side-effects.
- Hands free treatment, no fatigues to therapists.
- Patients can remain dressed up.
- 10-20 minutes long treatment times.





**Figure 1**  
Homato Applicator



**Figure 2**  
Homato Analogue Instrument







## Connection

Homato Applicator to the instrument.

# homato